



PORTABLE B-ONE

with carry case

INDUSTRIAL VISCOMETERS

The PORTABLE B-ONE is the ideal tool for your viscosity measurements in the workshop, production area, delivery area or outside. With one hour of battery life, it will optimize your control times and flexibility.



Charge indicator.

Autonomy more than one hour.

Easy handling.

Carrying case included.

PORTABLE B-ONE

SPECIFICATIONS

Activity domains



Food Industry



Cosmetics pharmaceuticals



Paint / ink / coatings



Chemical / petroleum products



Car industry



Building materials



Teaching



Type of instrument

Rotating springless viscometer with 7" Touch screen

Rotation speeds

Unlimited number of speeds between 0.3 and 250 rpm

Torque range

From 0.05 to 13 mNm

Accuracy

+/- 1 % of the full scale

Repeatability

+/- 0,2 %

Display

Viscosity – Speed – Torque – Time – Charge level indication
Choice of viscosity units: cP/Poises or mPa.s / Pa.s

Security and confidentiality

An «operator» function allows you to enter a username for you instrument. This user must then be identified using a 4-digit code. There is also a protected mode that locks your measurement conditions.

Language

French/English/Russian/Spanish

Compatible measuring system

MS ASTM, MS BV, MS VANE

Compatible temperature control

EVA LR-BV

Supply voltage

90-240 VAC 50/60 Hz

What benefits are there for you?

The Portable B-ONE gives you with an hour of measuring time. Set measuring times for your thixotropic products. Get kinematic viscosity by inputting the density of your product.

Carry case

Included

Dimensions and weight

Head: Ø 85 mm
Height: 310 mm
Box: L265 x W125 x H65 mm
Weight: 2 kg

Available instrument :

Part Number Instrument

T600700

Désignation Instrument

PORTABLE B-ONE VISCOMETER IN CARRYING CASE

Tel : +33 (0)4 78 08 54 06
Fax : +33 (0)4 78 08 69 44
contact@lamyreology.com

www.lamyreology.com

11 A, rue des Aulnes
69410 Champagne au Mont d'Or
France

